

RIDX™ ASFV Ag Test Kit

African swine fever (ASF) is a lethal hemorrhagic fever of pigs with high mortality and morbidity rates. Currently, no treatments or vaccines are available to prevent ASF, so the control strategy mainly relies on the strict administration of sanitary measures and slaughtering of infected and exposed animals. Therefore, rapid and accurate diagnosis of the infections is essential for preventing, controlling, and exterminating of the disease in the fields.

The RIDX™ ASFV Ag Test Kit is a lateral flow chromatographic immunoassay for the detection of african swine fever virus (ASFV) antigen in porcine blood.

Principal: Chromatographic immunoassay for the qualitative test
Specimen: Porcine whole blood

Components

- ASFV Ag test device (10 tests)
- Dilution buffer (1 vial)
- Disposable capillary tube (10 ea)

Features

Sensitivity 83.33% (10/12) vs. PCR

Specificity 100% (435/435) vs. PCR


Limit of Detection 5×10^4 TCID₅₀/ml, 7.5 HAD₅₀/ml

No Cross-reactivity with other porcine infectious pathogens (CSFV, JEV, PCV2, PRRSV, PEDV, and TGEV)




Test Procedure


- 1 Collect the sample (porcine anti-coagulated whole blood) using a capillary tube into the line on the capillary tube.




- 2 Add 20ul of the sample into the sample hole.



- 3 Add 3 drops of the dilution buffer into the sample hole.



- 4 Read the result at 15 minutes.



CAT No.	Description	LIC No.	Type	Package size
LGM-PAG-11	RIDX ASFV Ag Test Kit	300-024	LFA*	10 Tests/Kit

* LFA: Lateral flow assay